IHDENT® ZirkonBridge ST

Einleitung

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IHDENT® ZirkonBridge ST

Introduction

IHDENT® ZirkonBridge ST dental zirconia blanks are machined from biocompatible yttria stabilized zirconia nano powder formed by cold isostatic pressing. At full density, the zirconia has high strength, excellent biocompatibility and optimum aesthetics. IHDENT® ZirkonBridge ST blanks can be machined for dental restorations by different CNC or manual milling systems. Applications include anterior and posterior, single tooth or bridge restorations.

Properties IHDENT® ZirkonBridge ST

Sintering Density (g/cm ³)	6.07
Bending Strength (MPa)	1200
Average Grain Size (µm)	0.4
Amount of monoclinic phase after accelerated aging (%)	<10
Translucency	27% (sample thickness 1mm)
Radioactivity (Bq/g)	<0.1

Tooth Preparation instructions

- A chamfer or shoulder (with 0.5 - 1.0mm) with a rounded inner angle is ideal.
- For occlusal surfaces or incisal areas, reduce by 0.8 - 1.5mm.
- Edge radius: R = 0.5
- Preparation angle: 3° - 6°
- For bridge preparations, check for parallelism and avoid undercuts.

Note: for full contour crowns, the wall thickness needs to be a minimum of 0.5mm

Directions for Use

- For CAD/CAM systems or when using manual methods, Ihdent Zirkonbridge ST dental zirconia blanks are milled to the required work-piece shape. The milling must add approximately 25% to the original size to compensate for shrinkage during the final sintering process.
- The exact shrinkage rate or expansion factor should be programmed into the CAM software, please check the exact individual shrinkage factor on the zirconia blank.
- For long bridges, please keep the supporting work-piece on the lingual or palatal side to prevent distortion of the bridge on cooling
- While using blanks with aluminum or plastic holders, please make sure to leave a residue of zirconia material of not less than 1.5 mm attached to the holder, which can be disposed of together with the holder after the milling process is complete.




Staining

(For detailed staining instructions, please refer to the Instructions of your color liquid).
For full contour crowns, please use IHDENT® Color Liquid with 16 shades. Painting by brush is recommended. In case that the work-piece is bathed in the stain, strong staining will result.
After staining, dry under an infrared lamp for 30 - 180 minutes according to the length and thickness of the work-piece.

Sintering

Max Heating rate: <25 °C per minute
Final temperature: 1480 °C
Holding time at final temperature: 2 hours
Cooling should be done without temperature control in the closed furnace. Never open the furnace before it has reached a temperature of less than 80 °C to avoid thermal shocks.

Sintering graph (Recommended)

Finishing

After milling, the work-piece can be further shaped with a handpiece and diamond burs, gently, at under 10000 rpm. The work-piece must always be wet ground to avoid overheating. It is recommended to sandblast with 50 micron aluminum oxide at 2 - 4 bar, to clean the work-piece and increase its surface density.

Available

IHDENT® ZirkonBridge ST	98.5 mm Blank with Standardgroove	14 mm REF 7690	18 mm REF 7691	22 mm REF 7692	25 mm REF 7693
-------------------------	-----------------------------------	----------------	----------------	----------------	----------------

IHDENT® ZirkonBridge ST	for AmmanGirrbach CNC, with Plastic frame	14 mm REF 7810	18 mm REF 7811	22 mm REF 7812	25 mm REF 7813
-------------------------	---	----------------	----------------	----------------	----------------

IHDENT® ZirkonBridge ST	95iger Blank - for Zirkonzahn®	14 mm REF 7820	18 mm REF 7821	22 mm REF 7822	25 mm REF 7823
-------------------------	--------------------------------	----------------	----------------	----------------	----------------

Legend

- Read instructions
- Expiration date
- For single use only
- LOT Charge number
- Keep in a dry place
- Store tightly closed keep

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Ihdent® and Zirkonzahn® are registered trademarks.