Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® Sutures

Sterile surgical sutures

Characteristics

IHDENT® Sutures are sterile, 75 cm long, individually packed sutures for surgical use. They are available in non-resorbable versions, silk and polyester. In addition a degradable type is available. IHDENT® Suture material (silk, polyester, fast degradable) is in generally used as approximation suture in dental applications (for intra-oral use). None of the types is intended to be used where there is a risk of having a biological effect or a risk to be wholly or mainly absorbed in the body, neither it can be used in cardiovascular surgery.

IHDENT® Sutures are easy to knot and handle, highly tenacious and tissue compatible. The polyester and silk sutures which are available in various thicknesses correspond to the requirements of the European Pharmacopoeia for sterile, non-resorbable sutures and the United States Pharmacopoeia for non-resorbable surgical sutures.

Use

IHDENT® silk sutures are made from specially braided natural silk threads. They are coloured black for improved visibility on the operation field. The tissue compatible impregnation makes the suture easy to use. The specially braided suture is particularly suitable for subcutaneous stitches, ligatures, purse-string ligatures and for closing wounds. The IHDENT® Sutures polyester sutures are dyed green to avoid confusion with silk (black). The hydrophobic coating results in above average tissue compatibility and is non-capillary. The excellent properties of polyester sutures make them suitable for ligatures, purse-string ligatures, vessel occlusion and for closing wounds.

Tissue compatibility

The sutures can elicit an initial inflammatory reaction in tissue, which is followed by a gradual encapsulation of the suture by fibrous connective tissue. Progressive degradation of the proteinaceous silk fibre may result in gradual loss of the suture's tensile strength over time.

Contra-indication

If allergies against silk threads are known, Polyester threads should be used instead.

Recommendations for use

The choice of non-resorbable IHDENT® Suture material to be used depends upon the condition of the patient, the experience of the doctor, the technique and the size and type of wound. Infected wounds should receive the necessary treatment before stitches are applied. Care should be taken that the sutures are not damaged during use. The properties of the suture may be impaired through damage to the structure e.g. by crushing with surgical instruments or tweezers.

Sterility

IHDENT® Suture materials (Silk, Polyester) are gamma sterilized. IHDENT® degradable suture material is sterilized with ethylene-oxide. As long as the package is not damaged or opened sterility is guaranteed. Damaged packs should be disposed of. Resterilization can alter the properties of the material.

Storage

The sterile sutures should not be stored above 25 °C and should be protected from moisture and direct sunlight. IHDENT® Sutures should be used prior to the expiration date. Keep out of reach of children.

Package sizes

Silk (12 in 1 box)

Wire	needle shape	needle length	REF
3.0	1/2 ▼	14	SK23014B4
3.0	1/2	16	SK23016B4
3.0	1/2	21	SK23021B4
3.0	3/8	16	SK23016F4
3.0	3/8	18	SK23018F4
3.0	3/8	21	SK23021F4
3.0	3/8	26	SK23026F4
4.0	1/2	14	SK24014B4
4.0	1/2	16	SK24016B4
4.0	1/2	21	SK24021B4
4.0	3/8	16	SK24016F4
4.0	3/8	18	SK24018F4
4.0	3/8	21	SK24021F4
4.0	3/8	26	SK24026F4

Polyester (12 in 1 box)

4.0	1/2 ▼	14	PB24014B4
4.0	1/2	16	PB24016B4
4.0	1/2	21	PB24021B4
4.0	3/8	16	PB24016F4
4.0	3/8	18	PB24018F4
4.0	3/8	21	PB24021F4
5.0	1/2	14	PB25014B4
5.0	1/2	16	PB25016B4
5.0	1/2	21	PB25021B4
5.0	3/8	16	PB25016F4
5.0	3/8	18	PB25018F4
5.0	3/8	21	PB25021F4

Fast degradable (12 in 1 box)

4.0	3/8 ▼	16	PGAR21640F4
4.0	3/8	18	PGAR21840F4
5.0	3/8	16	PGAR21650F4
5.0	3/8	18	PGAR21850F4

Legends

Read instructions

Expiration date

LOT Charge number

Store protected against the sun

Keep in a dry place

Temperatures up to 25 °C

Gamma-sterilized

This product is only for sale to dentists, oral surgeons and maxillo-facial surgeons or on their behalf.

Sterilization by ethyleneoxide

Attention, see instructions for use

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IHDENT® is a registered trademark.
Dr. Ihde Dental AG
Dr. Ihde Dental GmbH
onewaybiomed GmbH