Deutsch:

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English:

Instructions for use – Please read carefully

IHDE[®] ZirkonBridge

Introduction

IHDE[®] ZirkonBridge ST dental zirconia blanks are machined from biocompatible yttria stabilized zirconia nano powder formed by cold isostatic pressing. At full density, the zirconia has high strength, excellent biocompatibility and optimum aesthetics. IHDE[®] ZirkonBridge ST blanks can be machined for dental restorations by different CNC or manual milling systems. Applications include anterior and posterior, single tooth or bridge restorations.

Contents

ZrO₂+HfO₂: 94 - 95 %
 Y₂O₃: 4.54 - 5.5 %
 Al₂O₃: 0.24 - 0.5 %
 Other oxides: <0.5 %

Properties IHDE[®] ZirkonBridge ST

Sintering Density (g/cm³): 6.07
 Bending Strength (MPa): 1200
 Average Grain Size (µm): 0.4
 Amount of monoclinic phase after accelerated aging (%): <10
 Translucency: 43 % (sample thickness 1 mm)
 Radioactivity (Bq/g): <0.1

Tooth Preparation instructions

1. A chamfer or shoulder (with 0.5 - 1.0mm) with a rounded inner angle is ideal.
 2. For occlusal surfaces or incisal areas, reduce by 0.8 - 1.5 mm.
 3. Edge radius: R = 0.5
 4. Preparation angle: 3° - 6°
 5. For bridge preparations, check for parallelism and avoid undercuts.
- Note:** for full contour crowns, the wall thickness needs to be a minimum of 0.5 mm.

Directions for Use

1. For CAD/CAM systems or when using manual methods, IHDE[®] ZirkonBridge ST dental zirconia blanks are milled to the required work-piece shape. The milling must add approximately 25 % to the original size to compensate for shrinkage during the final sintering process.
2. The exact shrinkage rate or expansion factor should be programmed into the CAM software, please check the exact individual shrinkage factor on the zirconia blank.
3. For long bridges, please keep the supporting work-piece on the lingual or palatal side to prevent distortion of the bridge on cooling
4. While using blanks with aluminum or plastic holders, please make sure to leave a residue of zirconia material of not less than 1.5 mm attached to the holder, which can be disposed of together with the holder after the milling process is complete.




Staining

The staining process depends very much on the kind of color. For detailed staining instructions, please refer to the instructions of your color liquid. Painting IHDE[®] ZirkonBridge ST Blanks by brush is recommended. In case that the work-piece is bathed in the stain, strong staining will result. After staining, dry under an infrared lamp for 30 - 180 minutes according to the length and thickness of the work-piece.

Sintering

Max Heating rate: < 25 °C per minute
 Final temperature: 1480 °C
 Holding time at final temperature: 2 hours
 Cooling should be done without temperature control in the closed furnace. Never open the furnace before it has reached a temperature of less than 80 °C to avoid thermal shocks.

Sintering graph (Recommended)

Finishing

After milling, the work-piece can be further shaped with a handpiece and diamond burs, gently, at under 10000 rpm. The work-piece must always be wet ground to avoid overheating. It is recommended to sandblast with 50 micron aluminum oxide at 2 - 4 bar, to clean the work-piece and increase its surface density.

Available:

IHDE[®] ZirkonBridge ST

98.5 mm Blank with Standardgroove
 12 mm REF 7689
 14 mm REF 7690
 18 mm REF 7691
 22 mm REF 7692
 25 mm REF 7693

IHDE[®] ZirkonBridge ST

for AmmanGirrbach CNC, with Plastic frame
 14 mm REF 7810
 18 mm REF 7811
 22 mm REF 7812
 25 mm REF 7813

IHDE[®] ZirkonBridge ST

95iger Blank - for Zirkonzahn[®]
 14 mm REF 7820
 18 mm REF 7821
 22 mm REF 7822
 25 mm REF 7823

Store away from children. For dental-practice use only.

Legend

-  Read instructions
-  Expiration date
-  For single use only
-  LOT Charge number
-  Keep in a dry place
-  Store tightly closed keep
-  Manufacturer

CE 0029

9999-0002-03

Revision issued 2015-08

Dr. Ihde Dental AG
 Dorfplatz 11
 CH-8737 Gommiswald/SG
 Tel. +41 (0)55 293 23 23
 Fax +41 (0)55 293 23 00
 contact@implant.com
 www.implant.com

Distribution by

Dr. Ihde Dental GmbH
 Erfurter Str. 19
 D-85386 Eching/München
 Tel. +49 (0)89 319 761-0
 Fax +49 (0)89 319 761-33
 info@ihde-dental.de
 www.ihde-dental.de

IHDE[®] and Zirkonzahn[®] are registered trademarks.

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